

COMPACT LAMELLAR FTIR SPECTROMETER

Far-IR THz Mid-IR Near-IR VIS



*light beam spectra
OPO OPA DFG
autocorrelation*

lamellar mirror technology
permanently aligned

The spectrometer L-FTS is ideal for continuous spectral monitoring of a collimated light beam, simultaneously from the visible to the far infrared in a single display, requiring no adjustments of its dispersion-free internal beamsplitter-optics-detection system, with a minimal footprint on the optical table.

The L-FTS features a novel permanently aligned lamellar mirror technology, and includes a mini-notebook for on-line Fourier transformation at up to 2 spectra/s as well as for averaging, ratioing and storing spectra. Internal sample holders enable transmission and reflection measurements with a user-supplied beam.

A user-supplied detector can be attached to the spectrometer. A second optical port allows external access to the interferometer for a user-supplied optical system for, e.g., double-pulse generation or autocorrelation scans to 7 ps.

- * ultra-wide coverage..... 0.4 – 300 μm
- * ultra-compact..... 15x12x7 cm^3
- * spectral resolution (nom.)..... 6.3 cm^{-1}
 - * with option R..... 3.4 cm^{-1}
- * clear aperture..... 17 mm dia.
- * intensity limit..... 2 W/cm^2
- * built-in window-less[‡] thermal detector:
 - * power limit..... 10 mW
 - * NEP (1Hz)..... 30 nW
 - * flat response..... 1 – 70 μm^{\ddagger}
 - * reduced response..... 0.4 – 1 μm
 - * reduced response..... 70 – 300 μm^{\ddagger}


[‡]A protective KBr window is available at no cost, however, it blocks wavelengths $>36 \mu\text{m}$.

For ordering write or call

LASNIX
Sonnenweg 32
82335 Berg
Germany
www.lasnix.com

Phone: +49 (0) 8151 953036
Fax: +49 (0) 91130844 88883
Email: info@lasnix.com

LASNiX LAMELLAR FTIR SPECTROMETER L-FTS




design principle

Au-coated lamellar mirrors (red) move through a second set of Au-coated mirrors (black), resulting in varying phase retardation: a wavefront-division interferometer
+ no interferometer arms, permanently aligned, drift-free
+ no dielectric involved, hence broadband VIS-IR-THz operation



measurement example

HeNe and CO₂ laser beams superimposed

interferogram



screenshot of spectrum pane



application example I

thermal IR from ca. 400°C wire
direct (black): Planck curve with air
absorption dips from H₂O and CO₂

same through plastic foil (red): extra
fingerprint absorption dips, and F.P.
thickness fringes

application example II

tuning supercontinuum
spectrum from nonlinear
fiber pumped by Er fiber
laser @ 1.55 μm

